



APPLANATION TONOMETER AT 900 & AT 900 D

Go for the Gold Standard

Tradition and Innovation – Since 1858, visionary thinking and a fascination with technology have guided us to develop innovative products of outstanding reliability: Anticipating trends to improve the quality of life.



Goldmann applanation tonometer

AT 900 and AT 900 D*

The gold standard in IOP measurement

Because it is fast and reproducible, Goldmann applanation tonometry has been the gold standard of tonometry for many decades until today. It is based on the Goldmann principle – proven during more than 50 years of usage with millions of patients.

Furthermore, thanks to the high precision of manufacturing at Haag-Streit the Goldmann applanation tonometer offers reliable performance during many years without the need for manufacturer maintenance.

AT 900 is the original version of the Goldmann tonometer and works 100% mechanically. The digital Goldmann tonometer AT 900 D* adds the advantages of state of the art digital technology.

Wireless transfer of the IOP readings and left/right eye position allows direct integration of the tonometer into an existing patient information system thus facilitating an ergonomic workflow.

*AT 900 D is not available in the US



TONOSAFE – THE DISPOSABLE TONOMETER PRISMS

Convenient – Safe – Accurate

Tonosafe avoids the need for timely and laborious disinfection without compromising the original Goldmann quality. It has been developed for single use only and eliminates therefore reliably the risk of cross infection even in the most difficult situations. Calibrated to Goldmann standards, Tonosafe delivers accurate results with all Goldmann and Perkins contact tonometers.

GOLDMANN TONOMETER PRISMS

Excellent optical quality

Goldmann measuring prisms distinguish themselves by their optical and mechanical performance which allows for fast and reliable measurements in daily practice. They are largely handmade by experienced workers and checked for low tolerances in weight, optical quality, durability and for consistent measurement results.






TONOSAFE



GOLDMANN TONOMETER PRISMS

Slit lamp model	Mod. R	Mod. T	Mod. BQ
BD 900	■		
BM 900	■	■	
BP 900		■	■
BQ 900		■	■



HAAG-STREIT AG

Gartenstadtstrasse 10
 3098 Koeniz
 Switzerland
 Phone +41 31 978 01 11
 Fax +41 31 978 02 82
 info@haag-streit.com
 www.haag-streit.com

Members of HAAG-STREIT Group

HAAG-STREIT Holding AG www.haag-streit-holding.com	HAAG-STREIT Deutschland GmbH www.haag-streit.de
HAAG-STREIT AG www.haag-streit.com	IPRO GmbH www.ipro.com
SPECTROS AG www.spectros.ch	CLEMENT CLARKE Ltd. www.clement-clarke.com
HAAG-STREIT Medtech AG www.haag-streit-medtech.com	HAAG-STREIT UK www.haag-streit-uk.com
HAAG-STREIT France EURL www.haag-streit.fr	John Weiss Ltd. www.johnweiss.com
HAAG-STREIT Far East www.haag-streit-fareast.com	HAAG-STREIT USA www.haag-streit-usa.com
HAAG-STREIT Surgical GmbH www.haag-streit-surgical.com	Reliance Medical Inc. www.haag-streit-usa.com
Möller-Wedel GmbH & Co KG www.moeller-wedel.com	Asetronics AG www.asetronics.ch
Möller-Wedel Optical GmbH www.moeller-wedel-optical.com	ComLab AG www.comlab.ch