





Tonoclear

Get a clear view of the mires during routine examinations, with our Goldmann-style, single-use applanation prism.



Tonoclear mires image



Competitor mires

Features

- **Superior image quality** – Tonoclear prisms are designed to be less affected by the intensity of light or internal reflections, enabling you to gain a clear image of the mires for improved diagnostic accuracy and speed of measurement.
- **Sterility** – Prisms are single-use and come ethylene oxide sterilized, to reduce the risk of cross-contamination between patients.
- **Pre-assembled for ease-of-use** – Individually packed and pre-assembled, making them quicker and easier to install, whilst further reducing the risk of cross-contamination.
- **Compatibility** – Works with most Goldmann applanation tonometers.



Single-use tonometer prisms for Goldmann applanation tonometry

Specification

Part number	2414-P-7045
Prism type	Single-use
Contents	100 x Tonoclear single-use prisms. Pre-assembled and individually packed.
Sterilization mode	Ethylene Oxide
Sterilization validity	3 years
Compatibility	Analogue Keeler Applanation Tonometers (KATs) – including R-Type, T-Type and Z-Type Digital Analogue Keeler Applanation Tonometers (D-KATs) – including R-Type and T-Type Digital Z-Type Keeler Applanation Tonometer (Z-KAT) Non-Keeler Goldmann Applanation Tonometers
Material	PMMA
Weight	1.65g ± 0.05
Prism diameter	6.90mm



Tonoclear prisms come pre-assembled in individual sterile blister packs, to reduce the risk of cross-contamination and make them quick and easy to install.



0325 REV A

Any questions? Please contact your local Keeler representative.